



## Contents

<b>1</b>	<b>Introduction</b>	<b>1</b>
<b>2</b>	<b>Prerequisites</b>	<b>1</b>
2.1	Register ZEN Functionality . . . . .	2
2.2	Starting MATLAB and ZEN . . . . .	3
<b>3</b>	<b>ZEN - Loading an CZI image</b>	<b>4</b>
<b>4</b>	<b>ZEN - Dialog to send the CZI to MATLAB via OAD</b>	<b>5</b>
<b>5</b>	<b>ZEN-MATLAB at work</b>	<b>7</b>
5.1	Importing CZI files into MATLAB . . . . .	9
5.2	Getting the CZI Metadata into MATLAB . . . . .	9
<b>6</b>	<b>Appendix: MATLAB M-Files</b>	<b>10</b>
6.1	ReadImage6D.m . . . . .	10
6.2	GetOMEDData.m . . . . .	11
<b>7</b>	<b>Disclaimer</b>	<b>13</b>

## 1 Introduction

This application note will explain how create a workflow using the ZEN-MATLAB connection. It illustrates the possibility to run MATLAB code from within a ZEN OAD macro.

## 2 Prerequisites

To be able to use the functionality described herein, you need

- **ZEN Blue 2012 or later with Macro Environment module**
- **MATLAB version R2013 - Standard**

To profit fully from the following example, you should have working experience with both applications. You should be familiar with the user interface in MATLAB, specifically with the Command Window and its use, and with the creation and use of scripts.



## 2.1 Register ZEN Functionality

To be able to use ZEN services in a .COM environment, provided by MATLAB, the ZEN functionality must be made registered as follows using a BAT-File, which **must** be executed as administrator.

```
1  echo off
2
3  pushd "C:\Windows\Microsoft.NET\Framework64\v4.0.30319"
4
5  SET dll-1="C:\Program Files\Carl Zeiss\ZEN 2\ZEN 2 (blue edition)\Zeiss.Micro.Scripting.dll"
6  regasm /u /codebase /tlb %dll-1%
7  regasm /codebase /tlb %dll-1%
8
9  SET dll-2="C:\Program Files\Carl Zeiss\ZEN 2\ZEN 2 (blue edition)\Zeiss.Micro.LM.Scripting.dll"
10 regasm /u /codebase /tlb %dll-2%
11 regasm /codebase /tlb %dll-2%
12 popd
13 pause
```

The final result of this set of commands is that the .NET classes and structures within **Zeiss.Micro.Scripting.dll** and **Zeiss.Micro.LM.Scripting.dll** are made known and accessible within .COM (**Component Object Model**) environment as well.

Remark: Please note that you must edit the BAT-file and adapt the path location, where the DLLs can be found to your needs.

## 2.2 Starting MATLAB and ZEN

This 1st thing one must do is to start MATLAB and check, if the automation services are activated.




Figure 1: MATLAB - Enable Automation Server

Now you can start the ZEN Blue software and create the following OAD Python script. Import is to check if the correct MATLAB Application ID was used.

### 3 ZEN - Loading an CZI image

Just load a CZI image of your choice into ZEN. For this example the selected image has the following dimensions:.

- Timepoints:  $T = 1$
- Z-Planes :  $Z = 22$
- Channels :  $CH = 3$
- Dimension Oder: XYZCT

Inside the ZEN Blue software this should similar tom this:




Figure 2: ZEN - The example image is loaded and displayed into ZEN

## 4 ZEN - Dialog to send the CZI to MATLAB via OAD

The main idea behind controlling ZEN from within MATLAB is just to make MATLAB aware of the ZEN Scripting DLLs via the .COM interface. This way you can use all the commands available inside the OAD python scripts directly within MATLAB and any M-file and vice versa.

This means you can either:

- Use OAD code directly within a M-File or from within the MATLAB command window
- Use MATLAB code directly within an OAD macro within ZEN

The following OAD script will show a very simple dialog that allows the user to select the dimension order and will send the active image to MATLAB using the second possibility mentioned above.

**Important Remark: MATLAB must be already running and the Automation Server must be activated.**

```
1  # Author: Sebastian Rhode
2  # Date: 26.03.2015
3  # Version: 1.2
4
5  # Send the current open image over to MATLAB using ReadImage6D.m --> has to be on the MATLAB path !!!
6  from System.Runtime.InteropServices import Marshal
7  from System.IO import Path, File, Directory, FileInfo
8  import os
9
10 CZIfiles_short = []
11 CZIdict = {}
12 # get all open documents
13 opendocs = Zen.Application.Documents
14 for doc in opendocs:
15     image = Zen.Application.Documents.GetByName(doc.Name)
16     if image.FileName.EndsWith('.czi'):
17         # get the filename of the current document only when it ends with '.czi'
18         CZIfiles_short.append(Path.GetFileName(image.FileName))
19         CZIdict[Path.GetFileName(image.FileName)] = image.FileName
20
21 ## Activate GUI
22 wd = ZenWindow()
23 wd.Initialize('Sent CZI to MATLAB',470,200,True,True)
24 ## add components to dialog
25 wd.AddLabel('Sends selected CZI image to MATLAB.', '0', '0')
26 wd.AddLabel('Uses MATLAB wrapper (bfoopen) for BioFormats.', '1', '0')
27 wd.AddDropDown('czi', 'Select CZI Image Data', CZIfiles_short, 0, '2','0')
28 ## show the window
29 result=wd.Show()
30 ## check, if Cancel button was clicked
31 if result.HasCanceled == True:
32     sys.exit('Macro aborted with Cancel!')
33
34 ## get the input values and store them
35 cziiname = result.GetValue('czi')
36
37 try:
38     MATLAB = Marshal.GetActiveObject('MATLAB.Application.8.5')
39     print 'ZEN-MATLAB bridge is OK.'
40     MLOK = True
41 except:
42     print 'MATLAB not running'
43
44 if MLOK == True:
```

```
45
46     ## get current active document
47     CZIfolder = os.path.dirname(CZIdict[cziname])
48     print 'Transfer: ', CZIdict[cziname]
49     ## create MATLAB variables
50     MATLAB.execute("filename = '"+CZIdict[cziname]+"")
51     MATLAB.execute("CZIimage = ReadImage6D(filename);")
52
53     print 'Done'
```

The OAD macro will open the following dialog:




Figure 3: ZEN - Display the dialog to send the image to MATLAB


The DropDown box contains all currently open CZI image documents. All the user has to do is to select the image to be send.the correct dimension order and press **OK**. The dialog then calls a number of MATLAB commands and used BioFormats to read the CZI into an 6D image array and the corresponding metadata into MATLAB

## 5 ZEN-MATLAB at work

The OAD macro behind this little dialog will execute just three lines of MATLAB code:


- filename = CZIfilename
- dimorder = DimOrder
- CZIimage = ReadImage6D(filename, dimorder)

Inside MATLAB the result will look like this:



As a result one will get a cell in MATLAB that contains two elements.

- a 6D matrix with the dimension [series, T, Z, C, X, X] depending on the dimension order
- a structure array containing metadata
  - **SeriesCount** = number of image series
  - **SizeC** = channels
  - **SizeT** = time points
  - **SizeZ** = focal planes
  - **SizeX** = pixels in X
  - **SizeY** = pixels in Y
  - **ScaleX** = scaling X in micron
  - **ScaleY** = scaling Y in micron
  - **ScaleZ** = spacing between two focal planes in micron
  - **DimOrder** = dimension order



The screenshot shows two MATLAB windows. The top window, titled 'Variables - CZImage{1, 2}', displays the fields of a structure array. The bottom window, titled 'Workspace', shows the variables in the current workspace.

Field	Value
SeriesCount	1
SizeC	3
SizeT	1
SizeZ	22
SizeX	425
SizeY	425
ScaleX	0.1062
ScaleY	0.1062
ScaleZ	0.3000
DimOrder	'XYCZT'

Name	Value
CZImage	1x2 cell
dimorder	'XYCZT'
filename	'C:\Users\M1SRH\Do...





## 5.1 Importing CZI files into MATLAB

ZEN Blue used the CZI data format to store the image data. And due to its open concept it is very well supported by the BioFormats library. Instructions on the usage of the MATLAB wrapper for BioFormats can be found here:

<http://www.openmicroscopy.org/site/support/bio-formats5/developers/matlab-dev.html>

There is more than one way to import a CZI file, so the choice depends on the nature of the data. For this example we use a pretty generic approach to read the CZI dataset. It is important to point out, the way the data were acquired has an influence on the way one must import them into MATLAB. Especially the order of the dimension CZT (Channels - Z -Planes - Time Points) can be important. The script used to read the CZI data can be found in section **6.1**

## 5.2 Getting the CZI Metadata into MATLAB

A very important topic are the CZI Metadata. Again we can rely on the BioFormats library and just use the already existing functionality to get all the information we need. For the example we use the following script shown in **6.2**.

## 6 Appendix: MATLAB M-Files

### 6.1 ReadImage6D.m

This is the M-file for reading then CZI data using the MATLAB toolbox for the BioFormats library.

```

1  % File: ReadImage6D.m
2  % Author: Sebastian Rhode
3  % Date: 19.09.2016
4  % Version: 1.2
5
6  % Read CZI image data into image6d array
7
8  function out = ReadImage6D(filename)
9
10 % Get OME Meta-Information
11 MetaData = GetOMEData(filename);
12
13 % The main inconvenience of the bfopen.m function is that it loads all the content of an image regardless of its size.
14 % Initialize BioFormats Reader
15 reader = bfGetReader(filename);
16
17 % add progress bar
18 h = waitbar(0,'Processing Data ...');
19 totalframes = MetaData.SeriesCount * MetaData.SizeC * MetaData.SizeZ * MetaData.SizeT;
20 framecounter = 0;
21
22 % Preallocate array with size (Series, SizeC, SizeZ, SizeT, SizeX, SizeY)
23 image6d = zeros(MetaData.SeriesCount, MetaData.SizeC, MetaData.SizeZ, MetaData.SizeT, MetaData.SizeX, MetaData.SizeY);
24
25 for series = 1: MetaData.SeriesCount
26
27     % set reader to current series
28     reader.setSeries(series-1);
29     for timepoint = 1: MetaData.SizeT
30         for zplane = 1: MetaData.SizeZ
31             for channel = 1: MetaData.SizeC
32
33                 framecounter = framecounter + 1;
34                 % update waitbar
35                 wstr = {'Reading Images: ', num2str(framecounter), ' of ', num2str(totalframes),'Frames' };
36                 waitbar(framecounter / totalframes, h, strjoin(wstr))
37
38                 % get linear index of the plane (1-based)
39                 iplane = loci.formats.FormatTools.getIndex(reader, zplane - 1, channel - 1, timepoint -1) +1;
40                 % get frame for current series
41                 image6d(series, timepoint, zplane, channel, :, :) = bfGetPlane(reader, iplane);
42
43             end
44         end
45     end
46 end
47
48 % close waitbar
49 close(h)
50
51 % close BioFormats Reader
52 reader.close();
53
54 % store image data and meta information in cell array
55 out = {};
56 % store the actual image data as 6d array
57 out{1} = image6d;
58 % store the image metainformation
59 out{2} = MetaData;
60

```

## 6.2 GetOMEData.m

This M-File reads all the relevant MetaData from the image. Please feel free to adapt it to your needs.

```

1 % File: GetOMEData.m
2 % Author: Sebastian Rhode
3 % Date: 13.09.2015
4 % Version: 1.3
5
6 function OMEData = GetOMEData(filename)
7 %
8 % Get OME Meta Information using BioFormats Library 5.1.10
9
10 % To access the file reader without loading all the data, use the low-level bfGetReader.m function:
11 reader = bfGetReader(filename);
12
13 % You can then access the OME metadata using the getMetadataStore() method:
14 omeMeta = reader.getMetadataStore();
15
16 % get ImageCount --> currently only reading one image is supported
17 imagecount = omeMeta.getImageCount();
18 % create empty cell array to store the image IDs
19 imageIDs_str = cell(1, imagecount);
20 imageIDs = cell(1, imagecount);
21
22 % try to get all the imageIDs as strings and numbers (zero-based)
23 try
24     for id = 1:imagecount
25         imageIDs_str{id} = omeMeta.getImageID(id-1);
26         imageIDs{id} = id-1;
27     end
28     % store the OMEData
29     OMEData.ImageIDs = imageIDs;
30     OMEData.ImageIDstrings = imageIDs_str;
31 catch
32     OMEData.ImageIDs = 'na';
33     OMEData.ImageIDstrings = 'na';
34     msg = 'No suitable ImageIDs found.';
35     warning(msg);
36 end
37
38 % use default imageID to read other metadata
39 imageID = imageIDs{1};
40
41 % get the actual metadata and store them in a structured array
42 [pathstr,name,ext] = fileparts(filename);
43 OMEData.FilePath = pathstr;
44 OMEData.Filename = strcat(name, ext);
45
46 % Get dimension order
47 OMEData.DimOrder = char(omeMeta.getPixelsDimensionOrder(imageID).getValue());
48
49 % Number of series inside the complete data set
50 OMEData.SeriesCount = reader.getSeriesCount();
51
52 % Dimension Sizes C - T - Z - X - Y
53 OMEData.SizeC = omeMeta.getPixelsSizeC(imageID).getValue();
54 OMEData.SizeT = omeMeta.getPixelsSizeT(imageID).getValue();
55 OMEData.SizeZ = omeMeta.getPixelsSizeZ(imageID).getValue();
56 OMEData.SizeX = omeMeta.getPixelsSizeX(imageID).getValue();
57 OMEData.SizeY = omeMeta.getPixelsSizeY(imageID).getValue();
58
59 % Scaling XYZ
60 try
61     OMEData.ScaleX = round(double(omeMeta.getPixelsPhysicalSizeX(imageID).value()),3); % in micron
62 catch
63     msg = 'Problem getting X-Scaling. Use Default = 1';
64     warning(msg);
65     OMEData.ScaleX = 1;
66 end
67
68 try
69     OMEData.ScaleY = round(double(omeMeta.getPixelsPhysicalSizeY(imageID).value()),3); % in micron
70 catch
71     msg = 'Problem getting Y-Scaling. Use Default = 1';
72     warning(msg);
73     OMEData.ScaleY = 1;
74 end
75
76
77 try
78     OMEData.ScaleZ = round(double(omeMeta.getPixelsPhysicalSizeZ(imageID).value()),3); % in micron

```

## Application Note - Control ZEN from MATLAB

```

79 catch
80     % in case of only a single z-plane set to 1 micron ...
81     msg = 'Problem getting Z-Scaling. Use Default = 1';
82     warning(msg);
83     OMEData.ScaleZ = 1;
84 end
85
86 % read relevant objective information from metadata
87 try
88     % get the correct objective ID (the objective that was used to acquire the image)
89     tmp = char(omeMeta.getInstrumentID(imageID));
90     OMEData.InstrumentID = str2double(tmp(end));
91     tmp = char(omeMeta.getObjectiveSettingsID(OMEData.InstrumentID));
92     objID = str2double(tmp(end));
93     % error handling --> sometime only one objective is there with ID > 0
94     numobj = omeMeta.getObjectiveCount(OMEData.InstrumentID);
95     if numobj == 1
96         objID = 0;
97     end
98
99     OMEData.ObjID = objID;
100 catch
101     msg = 'No suitable instrument and objective ID found.';
102     warning(msg);
103 end
104
105 try
106     % get objective immersion
107     OMEData.ObjImm = char(omeMeta.getObjectiveImmersion(OMEData.InstrumentID, OMEData.ObjID).getValue());
108 catch
109     msg = 'Problem getting immersion type.';
110     warning(msg);
111     OMEData.ObjImm = 'na';
112 end
113
114 try
115     % get objective lens NA
116     OMEData.ObjNA = round(omeMeta.getObjectiveLensNA(OMEData.InstrumentID, OMEData.ObjID).doubleValue(),2);
117 catch
118     msg = 'Problem getting objective NA.';
119     warning(msg);
120     OMEData.ObjNA = 'na';
121 end
122
123 try
124     % get objective magnification
125     OMEData.ObjMag = round(omeMeta.getObjectiveNominalMagnification(OMEData.InstrumentID, OMEData.ObjID).doubleValue(),2);
126 catch
127     msg = 'Problem getting objective magnification.';
128     warning(msg);
129     OMEData.ObjMag = 'na';
130 end
131
132 try
133     % get objective model
134     OMEData.ObjModel = char(omeMeta.getObjectiveModel(OMEData.InstrumentID, OMEData.ObjID));
135 catch
136     msg = 'Problem getting objective model.';
137     warning(msg);
138     OMEData.ObjModel = 'na';
139 end
140
141 % get excitation and emission wavelengths for all channels
142 for c = 1:OMEData.SizeC
143     try
144         OMEData.WLEx{c} = round(omeMeta.getChannelExcitationWavelength(imageID, c-1).value().doubleValue());
145         OMEData.WLEm{c} = round(omeMeta.getChannelEmissionWavelength(imageID, c-1).value().doubleValue());
146     catch
147         msg = 'Problem getting excitation and emission wavelengths. Set to zero.';
148         warning(msg);
149         OMEData.WLEx{c} = 0;
150         OMEData.WLEm{c} = 0;
151     end
152
153     try
154         OMEData.Channels{c} = char(omeMeta.getChannelName(imageID, c-1));
155         OMEData.Dyes{c} = char(omeMeta.getChannelFluor(imageID, c-1));
156     catch
157         msg = 'No Metadata for current channel available.';
158         warning(msg);
159         OMEData.Channels{c} = 'na';
160         OMEData.Dyes{c} = 'na';
161     end
162 end
163
164 % close BioFormats Reader
165 reader.close()
166

```



## 7 Disclaimer

This is an application note free to use for everybody. Use it on your own risk.

Carl Zeiss Microscopy GmbH's ZEN software allows to connect to a the third party software MATLAB. Therefore Carl Zeiss Microscopy GmbH undertakes no warranty concerning MATLAB, makes no representation that MATLAB will work on your hardware and will not be liable for any damages caused by the use of this example. By running this example you agree to this disclaimer.